



HY30-D
Wrist Laser Treatment
Instrument



Features:

- Laser irradiation to three parts or acupoints of human simultaneously, it's an unprecedented function in the market.
- A laser stabilizer specially added to perform with more stability.
- 6 pieces national patent protected.
- Proper design to fit human physiological curve and ergonomics.
- High capacity lithium battery for long-term use.
- Adjustable treatment time/power, intelligent memory system.
- Two output modes: Pulse and Continues, can be switched manually or automatically.

HY30-D

Wrist Laser Treatment
Instrument



Application Range:

Hyperviscosity, hyperlipidemia, hypertension, diabetes, ischemic cardiovascular and cerebrovascular diseases. It's also has good curative on the Rhinitis, allergic rhinitis, acute rhinitis, chronic rhinitis, sinusitis and nasal polyps.

Specification:

Item	Description
Laser medium	GaAlAs Semiconductor
Laser wavelength	650nm
Laser output power	≤5mw
Terminate laser medium quantity	6pcs.
Terminate output laser control	3 grade for adjustment
Timing range	10 ~ 60min, six grades for adjustment
Environment temperature	5°C-40°C
Relative humidity	≤80%
Atmospheric pressure	86kpa-106kpa
Battery voltage	3.7 to 4.2V
Lithium ion battery output voltage	5.5V
Power source use	lithium ion battery with the capacity 1600m Ah.
Output frequency under Pulse mode	0.5sec
Instrument power consumption	< 1 W

